



## **BLOOD GLUCOSE: (CAPILLARY MONITORING, BEDSIDE)**

**Choose the correct answer.**

1. A Protective Sleeve must be used to cover the monitor for all patients in any type of isolation and for all neonates to prevent cross contamination.
  - a. True
  - b. False
  
2. The needle of the lancing device must be disposed of in the sharps container after each use.
  - a. True
  - b. False
  
3. When performing a patient test, if the Precision PCx fails to autostart, how long do you have to apply a second drop of blood after the first drop of blood?
  - a. 15 seconds
  - b. 30 seconds
  - c. 60 seconds
  - d. a second drop can not be applied
  
4. How often must quality controls be performed on the monitor?
  - a. every shift
  - b. every 24 hours
  - c. once a week
  - d. once a month
  
5. If a bottle of control solution is opened on 1-1-99, it will expire:
  - a. 30 days after it is opened
  - b. 90 days after it is opened
  - c. on the date indicated on the bottle
  
6. A blood specimen will be sent to the laboratory for comparison testing when critical values of:
  - a. >600mg/dl or <20 mg/dl for adults
  - b. 40 mg/dl or less for pediatric service line
  - c. >500 mg/dl or <40 mg/dl
  - d. a & b

## BLOOD GLUCOSE

7. For every patient that had a capillary blood glucose one must:
  - a. Order a GID
  - b. Call the lab
  - c. Send the GID to the lab as soon as possible
  - d. A & B

\*\*\*End\*\*\*